

Development of rapid *Legionella* viability qPCR detection assay to reduce risk of disease outbreaks associated with cooling towers

Dr. Nishanth Kuchi¹, Dr. Attiq Rehman^{1,2}, Matthew Poole¹, Nadine Godin¹, Michael Lawlor¹, Dr. Diane Botelho^{1,2} ¹Research and Productivity Council (RPC), 115A Harrisville Blvd, Moncton NB Canada E1H 3T3; ²Research and Productivity Council (RPC), 921 College Hill Rd, Fredericton NB Canada E3B 6Z9 Contact: nishanth.kuchi@rpc.ca; attiq.rehman@rpc.ca

INTRODUCTION

Legionellosis has been an increasing public health concern in New Brunswick. The severe disease condition is predominantly caused by the bacterium *Legionella pneumophila (LP)* which thrives in moist environments such as cooling towers. Traditional culture-based detection approaches are time-consuming, since results take 7–9 days. These approaches also lack the ability to detect "viable but non-culturable" *Legionella* bacteria that are often residually present after common disinfection treatments (e.g., chlorination or heat). Quantitative real time PCR (qPCR) based detection is a rapid molecular method that directly addresses these issues for accurate identification of *LP*. The main drawback of existing qPCR approaches is the lack of live-dead discrimination of the bacterium. We've developed a *Legionella* viability qPCR detection method which combines the rapidity of qPCR with the ability to discriminate between viable and non-viable cells (Figure 1). Our objective is to develop a high-throughput, rapid qPCR detection method for *LP*, *LP* serogroup 1 (LPSG1) and other *Legionella* species.

MATERIALS & METHODS

Water samples were spiked with known quantities (High – 10^3 to 10^4 MPN/mL; Medium – 10^2 to 10^3 MPN/mL; Low – 10^1 to 10^2 MPN/mL) of target (*LP*) and non-target (*E.coli*) viable and non-viable pathogens. The spikes were filter concentrated and subjected to DNA extraction followed by qPCR analysis¹ on the Bio-Rad Opus CFX-96 platform. The method contains a cell-permeable intercalating crosslinker that exclusively binds to DNA from membrane-compromised dead, dying, or lysed cells^{2,3} (Figure 1). Once bound, the molecule irreversibly crosslinks DNA, making it un-amplifiable in PCR reactions. As a result, only the signal exclusively from viable cells is detected. The data was analyzed using Bio-Rad CFX Maestro software (Figures 2a, b) and qPCR results were compared with traditional methods⁴. Relative quantification was performed using a DNA standard with known copies of target DNA.



RESULTS AND DISCUSSION

The qPCR protocol was highly specific to *Legionella* spp. (Table 1) and differentiated from non-target (*E.coli*) and negative controls (no amplification was observed). The results based on inclusivity and exclusivity show that the method can accurately distinguish between viable and non-viable cells (Figure 3).

Acceptable	Efficiency	Slope	R²
Range	(86-115%)	(-3.0 to -3.7)	(>0.970)
Observed range	89.57 to 100.07	-3.32 to -3.6	>0.995

Table 1: Evaluation of the qPCR Standard Curve

Spike	MPN±RSD/mL	vGE ±RSD/mL
High	> 2272.6	1791.4 ± 28.4
Medium	970.9 ± 23.3	116.5 ± 33.0
Low	114.1 ± 40.9	24.3 ± 113.5

Table 2: Comparison of Culture-Based Enumeration of LP via Most ProbableNumber (MPN/mL) and Viable Genomic Equivalent qPCR vGE/mL

rpc







Figure 3: Viable (treated with crosslinker (CL)) and total (without CL) gene copies in live and dead samples of a) *Legionella* spp.; b) *L. pneumophila*; and c) *L. pneumophila* Serogroup 1 in high (blue), medium (orange) and low (grey) spikes. GE/mL – Genomic equivalent/mL; vGE/mL – viable genomic equivalent/mL. Error bars represent RSD%.

The limit of detection and quantification of the assay was 20 copies/reaction. Based on the results, vGU/mL can be used as an accurate, rapid and informative alternative to the more commonly used culture-based methods.

CONCLUSIONS

- The viability qPCR method is a rapid and reliable alternative to traditional culture-based methods for the detection of Legionella spp., L. pneumophila and L. pneumophila Serogroup 1.
- This rapid detection method (~3-4 hours) can be used to assist public health/environmental professionals in performing risk assessment and allows for rapid response to mitigate outbreaks.

ACKNOWLEDGEMENTS

This study was supported by the Total Development Funds (TDF) from Regional Development Corporation (RDC). We would like to thank Denis Chenard (Healthy Environments Branch, Department of Environment and Local Government). Our colleagues from RPC, Dr. Tony Manning, Emma Matchet and RPC's Business Growth and Excellence Team including Josh Perry, Amy Brown, Scott Archibald and Sara Cockburn for their contributions to this project. We thank Dr. Nathan Feirer (Senior Research Scientist) and Dr. Carl Strayer (Technical Services Quality Specialist), Promega, Madison, WI, USA for their help during the standardization.

REFERENCES

1. Promega GoTaq Legionella vPCR kit protocol TM#725, Promega, Madison, WI, USA.

2. Toplitsch, D., Platzer, S., Zehner, R., Maitz, S., Mascher, F., & Kittinger, C. (2021). Comparison of updated methods for Legionella detection in environmental water samples. International Journal of Environmental Research and Public Health, 18(10), 5436.

3. Delgado-Viscogliosi, P., Solignac, L., & Delattre, J. M. (2009). Viability PCR, a culture-independent method for rapid and selective quantification of viable Legionella pneumophila cells in environmental water samples. Applied and environmental microbiology, 75(11), 3502-3512..

4. Lizana, X., López, A., Benito, S., Agustí, G., Ríos, M., Piqué, N., Marques, A.M., & Codony, F. (2017). Viability qPCR, a new tool for Legionella risk management. International Journal of Hygiene and Environmental Health, 220(8), 1318-1324.